

**PRODUCT LIST**  
**OCCLUTECH AB-SWEDEN**

***PDA OCCLUDER***

*D2   D3   D1   Device Length (mm)*  
*mm   mm   mm*

<i>Product No.</i>	<i>Delivery System(F)</i>				
<i>Standard Shank</i>					
<b>42PDA05</b>	<b>6</b>	<b>3.5</b>	<b>5</b>	<b>9</b>	<b>4.25</b>
<b>42PDA06</b>	<b>6</b>	<b>4</b>	<b>6</b>	<b>10</b>	<b>5.00</b>
<b>42PDA07</b>	<b>6</b>	<b>5</b>	<b>7</b>	<b>11</b>	<b>6.05</b>
<b>42PDA08</b>	<b>6</b>	<b>6</b>	<b>8</b>	<b>13</b>	<b>6.30</b>
<b>42PDA10</b>	<b>7</b>	<b>8</b>	<b>10</b>	<b>16</b>	<b>7.00</b>
<b>42PDA12</b>	<b>7</b>	<b>10</b>	<b>12</b>	<b>18</b>	<b>12.00</b>
<b>42PDA15</b>	<b>8</b>	<b>12</b>	<b>15</b>	<b>20</b>	<b>14.00</b>
<b>42PDA18</b>	<b>9</b>	<b>14</b>	<b>18</b>	<b>24</b>	<b>16.00</b>
<i>Long Shank</i>					
<b>43PDA05L</b>	<b>6</b>	<b>3.5</b>	<b>5</b>	<b>9</b>	<b>7.00</b>
<b>43PDA06L</b>	<b>6</b>	<b>4</b>	<b>6</b>	<b>10</b>	<b>7.50</b>
<b>43PDA07L</b>	<b>6</b>	<b>5</b>	<b>7</b>	<b>11</b>	<b>8.50</b>
<b>43PDA08L</b>	<b>6</b>	<b>6</b>	<b>8</b>	<b>13</b>	<b>9.00</b>
<b>43PDA10L</b>	<b>7</b>	<b>8</b>	<b>10</b>	<b>16</b>	<b>10.50</b>

*D1: Retention disk   D2: Aortic Side of the Shank   D3: Pulmonic side Of the Shank*

**OCCLUTECH ASD Occluder**

<i>Product No.</i>	<i>Waist mm</i>	<i>LA Disc mm</i>	<i>RA Disc mm</i>
<b>29ASD04</b>	<b>4</b>	<b>11</b>	<b>9</b>
<b>29ASD05</b>	<b>5</b>	<b>14</b>	<b>11</b>
<b>29ASD06</b>	<b>6</b>	<b>16.5</b>	<b>12.5</b>
<b>29ASD07</b>	<b>7.5</b>	<b>18</b>	<b>14</b>
<b>29ASD09</b>	<b>9</b>	<b>20.5</b>	<b>16.5</b>
<b>29ASD10</b>	<b>10.5</b>	<b>22</b>	<b>18</b>
<b>29ASD12</b>	<b>12</b>	<b>27</b>	<b>23</b>
<b>29ASD13</b>	<b>13.5</b>	<b>28.5</b>	<b>24.5</b>
<b>29ASD15</b>	<b>15</b>	<b>30</b>	<b>26</b>
<b>29ASD16</b>	<b>16.5</b>	<b>31.5</b>	<b>27.5</b>
<b>29ASD18</b>	<b>18</b>	<b>33</b>	<b>29</b>
<b>29ASD19</b>	<b>19.5</b>	<b>34.5</b>	<b>30.5</b>
<b>29ASD21</b>	<b>21</b>	<b>36</b>	<b>32</b>
<b>29ASD24</b>	<b>24</b>	<b>39</b>	<b>35</b>
<b>29ASD27</b>	<b>27</b>	<b>42</b>	<b>38</b>
<b>29ASD30</b>	<b>30</b>	<b>45</b>	<b>41</b>
<b>29ASD33</b>	<b>33</b>	<b>48</b>	<b>43</b>
<b>29ASD36</b>	<b>36</b>	<b>52</b>	<b>46</b>
<b>29ASD39</b>	<b>39</b>	<b>54</b>	<b>49</b>
<b>29ASD40</b>	<b>40</b>	<b>55</b>	<b>50</b>

*OCCLUTECH Mysd Occluder*

<b>Product No.</b>	<b>Height mm</b>	<b>D1 mm</b>	<b>D2 mm</b>
<b>71VSD04</b>	<b>7</b>	<b>4</b>	<b>10</b>
<b>71VSD06</b>	<b>7</b>	<b>6</b>	<b>13</b>
<b>71VSD08</b>	<b>7</b>	<b>8</b>	<b>15</b>
<b>71VSD10</b>	<b>7</b>	<b>10</b>	<b>17</b>
<b>71VSD12</b>	<b>7</b>	<b>12</b>	<b>19</b>
<b>71VSD14</b>	<b>7</b>	<b>14</b>	<b>22</b>
<b>71VSD16</b>	<b>7</b>	<b>16</b>	<b>24</b>
<b>71VSD18</b>	<b>7</b>	<b>18</b>	<b>26</b>
<b>71VSD20</b>	<b>7</b>	<b>20</b>	<b>28</b>

*D1-Diameter of the Waist*

*D2-Diameter of the discs*

*H-Height of the waist*

**OCCLUTECH PFO Occluder**

<b>Product No.</b>	<b>LA Disc mm</b>	<b>RA Disc mm</b>
<b>LA Single Layer</b>		
<b>18PFO255</b>	<b>23</b>	<b>25</b>
<b>LA Double Layer</b>		
<b>19PFO18D</b>	<b>16</b>	<b>18</b>
<b>19PFO25D</b>	<b>23</b>	<b>25</b>
<b>19PFO30D</b>	<b>27</b>	<b>30</b>
<b>19PFO35D</b>	<b>31</b>	<b>35</b>

**OCCLUTECH Sizing Balloon**

<b>Product No</b>	<b>Balloon Nominal mm</b>	<b>Max Fill Vol(mm)</b>	<b>Balloon Length mm</b>
<b>53SB254</b>	<b>25</b>	<b>35</b>	<b>27</b>
<b>53SB355</b>	<b>35</b>	<b>90</b>	<b>40</b>